

## RADCAL TOUCH SERIES COMPATIBILITY AND FEATURES TABLE:

|                                 | Accu-Gold+ |           | Rapid-Gold+ |           | Accu-Dose+ |           |
|---------------------------------|------------|-----------|-------------|-----------|------------|-----------|
|                                 | Touch      | Touch Pro | Touch       | Touch Pro | Touch      | Touch Pro |
| <b>Functional Compatibility</b> |            |           |             |           |            |           |
| Dose                            | •          | •         | •           | •         | •          | •         |
| Dose Rate                       | •          | •         | •           | •         | •          | •         |
| Exposure Time                   | •          | •         | •           | •         | •          | •         |
| Pulse Count                     | •          | •         | •           | •         | •          | •         |
| kV                              | •          | •         | •           | •         |            |           |
| HVL                             | •          | •         | •           | •         |            |           |
| Filtration                      | •          | •         | •           | •         |            |           |
| Current (mA)                    | •          | •         | •           | •         |            |           |
| Luminance and Illuminance       |            | •         |             |           |            |           |
| Scatter and Leakage             | •          | •         |             |           | •          | •         |
| <b>Sensor Compatibility</b>     |            |           |             |           |            |           |
| Accu-Gold Multi-Sensors (AGMS)  | •          | •         | •           | •         |            |           |
| Accu-Gold Ion Chambers (10x6)   | •          | •         |             |           | •          | •         |
| Accu-Gold Dose Diodes (DDX6)    | •          | •         | •           | •         | •          | •         |
| mA and mAs Sensore (90M)        | •          | •         | •           | •         |            |           |
| Light Sensor (AGLS)             |            | •         |             |           |            |           |
| <b>Computer Features</b>        |            |           |             |           |            |           |
| Accu-Gold Windows Application   |            | •         |             | •         |            | •         |
| Waveforms                       |            | •         |             | •         |            | •         |
| Excel Reporting                 |            | •         |             | •         |            | •         |
| iOS and Android Applications    |            | •         |             | •         |            | •         |
| USB Connectivity                |            | •         |             | •         |            | •         |
| Wireless Connectivity           |            | •         |             | •         |            | •         |