

DDX6-WL

Diagnostic range
(Low dose rates)



for use with: *Accu-Gold+ / Rapid-Gold+ / Accu-Dose+*

ACCU-GOLD+

RAPID-GOLD+

ACCU-DOSE+

SPECIFICATIONS / TECHNICAL DATA:

All specifications subject to change.

SENSORS:	DDX6-WL	
Min Rate:	70 nR/s	600 pGy/s
Max Rate:	640 mR/s	5.6 mGy/s
Min Dose:	410 nR	4 nGy
Max Dose:	1.7 kR	15 Gy
Accuracy:	±5% using X-rays @ 80kVp with 2.5 mm Al total filtration (IEC 61267 RQR-6)	
Energy dependence:	± 5% from 50 kVp to 120 kVp at 2.5 mm Al	
Filtration dependence:	+10% to -5% from 2.5 mm Al to 23 mm Al	