

# Appendix D

## **Declaration of Conformity**

*According to Test Report No. 97339-4 and 97339-5*

**Manufacturer's Name:** Radcal  
**Manufacturer's Address:** 426 West Duarte Road  
Monrovia, CA 91016  
U.S.A.

## **Declares that the Product**

**Product Name:** Accu-Gold Wi-Fi Module  
**Model Number(s):** AGNUGGET/AGNUGGET-JP  
**Product Accessories:** All

## **Conforms to the following Product Specification:**

**EMC:** EN61326-1 (2006) Group 1 Class A  
EN 61000-4-2 (1995) +A1 +A2  
EN 61000-4-3 (2002)  
EN 61000-4-4 (2004)  
EN 61000-4-6 (2009)

The product herewith complies with the RFI-emissions<sup>1</sup> requirements and immunity requirements and carries the CE marking accordingly.

E. MacIntosh  
Engineering  
10 Aug 2015

---

<sup>1</sup> This product is known to cause interference when used with the charger connected and powered. If interference occurs in a particular installation the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the instrument.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Use it without the charger.