



## Accu-Gold 2.0 - Download

### Accu-Gold 2.0 (Windows Software)

- Software Setup V2.75.3
- ReadMe - check here for installation notes and cor
- Sample Templates

Verify: AG2/Help/About:

Accu-Gold 2

Version 2.75 Build 3

✓ Accu-Gold is up to date.


## Accu-Gold Touch – Download

(Flash Drive setup or  
Zip download)

### Accu-Gold Touch Firmware

- AGT Update V1.60.0 (Flash Drive setup)
- AGT Update V1.60.0 (Alt - Zip download)
- Firmware upgrade instructions

Verify:





## AGNugget WiFi – Download

(Flash Drive setup or  
Zip download)

### Accu-Gold Nugget Firmware

- Nugget Update V1.25.0 (Flash Drive setup)
- Nugget Update V1.25.0 (Alt - Zip download)
- Firmware upgrade instructions

Verify: Using Nugget, after making measurement, check Sensor tab

Digitizer	
Sensor ID	40-0002
CPC-Date	6/11/2020
Application Version	2.75 Build 3
Connection Type	Wi-Fi
Firmware Version	01.00
Nugget Firmware Version	1.25.0.0
Nugget ID	Nugget_51-0012
Nugget Wi-Fi MAC	00:25:97:5E:10:AB
Nugget Wi-Fi Version	2.6.2,1.7.7

## AGTouch WiFi – Download

(Flash Drive setup or  
Zip download)

### Accu-Gold Touch WiFi Firmware

- AGTWiFi Update V1.25.0 (Flash Drive setup)
- AGTWiFi Update V1.25.0 (Alt - Zip download)
- Firmware upgrade instructions

Verify: Using Touch WiFi, after making measurement, check Sensor tab

Digitizer	
Sensor ID	48-1090
CPC-Date	10/8/2019
Application Version	2.75 Build 3
Connection Type	Wi-Fi
Firmware Version	1.51
Nugget Firmware Version	1.25.0.0
Nugget ID	Touch_55-1090
Nugget Wi-Fi MAC	00:25:97:0C:06:A0
Nugget Wi-Fi Version	2.6.2,1.7.7

Note that the version numbers for the Nugget and the Touch WiFi will be the same but the “NuggetID” will identify one or the other.