

	Accu-Gold+	Rapid-Gold+	Accu-Dose+
Functional Compatibility			
Dose	●	●	●
Dose Rate	●	●	●
Waveform	●	●	●
Exposure Time	●	●	●
kV	●	●	
HVL	●	●	
Filtration	●	●	
Current (mA)	●	●	
Luminance and Illuminance	●		
Scatter and Leakage	●		●
Sensor Compatibility			
Accu-Gold Multisensors (AGMS)	●	●	
Accu-Gold Ion Chambers (10X6)	●		●
Accu-Gold Dose Diodes (DDX6)	●	●	●
mA and mAs Sensors (90M)	●	●	
Light Sensor (AGLS)	●		
DAP Check Plus Sensor (DAPCKP+)	●		