

The Rapid-Gold+ is a highly capable diagnostic measurement solution offering support for Radcal's Solid-State sensors and advanced diagnostic functionality.

Systems Solutions

The Rapid-Gold+ provides a tailored solution to your individual diagnostic measurement needs when paired with:

- Radcal's broad selection of Gold Standard sensors
- user friendly display that facilitates recording and reporting of results quickly and easily
- carry case to transport your custom solution securely and conveniently

Sensor Selection

The Rapid-Gold+ supports Radcal's full line of Solid-state Dose Diodes, and Solid-state Multisensors featuring the most accurate and compact stacked sensor design available.

Software and Display Technology

The Rapid-Gold+ offers great flexibility in your choice of user interface. Accu-Gold and Accu-Gold Excel Windows software provides a rich user-friendly environment in which to record, view, and archive your measurements. Radcal offers a tablet computer to run this software or use your own Windows computer. Apps are also available for your Apple iOS and Android devices (Nugget wireless adapter required). Alternatively, upgrade your instrument to an Rapid-Gold+ Touch stand-alone system equipped with a built-in display.

Functionality

The Rapid-Gold+ provides a comprehensive set of parameters including Dose, Dose Rate, Waveform, Exposure time, kV, Filtration, HVL, and mA (optional).

Modalities

The Rapid-Gold+ is well suited for Radiography, Fluoroscopy, Mammography, CT, and Dental applications



Unparalleled Accuracy, Experience, Versatility



426 West Duarte Road Monrovia, CA 91016 USA T: (626) 357-7921 F: (626) 357-8863 www.radcal.com

RAPID-GOLD+ FEATURES, BENEFITS AND SENSORS:



Full Specifications & Further Information: www.radcal.com