

RADCAL TOUCH SERIES COMPATIBILITY AND FEATURES TABLE:

Model Number	Accu-Gold+		Rapid-Gold+		Accu-Dose+	
	T3 AGT3-AG	T3 Pro AGT3-P-AG	T3 AGT3-RG	T3 Pro AGT3-P-RG	T3 AGT3-AD	T3 Pro AGT3-P-AD
Functional Compatibility						
Dose	•	•	•	•	•	•
Dose Rate	•	•	•	•	•	•
Exposure Time	•	•	•	•	•	•
Pulse Count	•	•	•	•	•	•
kV	•	•	•	•		
HVL	•	•	•	•		
Filtration	•	•	•	•		
Current (mA)	•	•	•	•		
Luminance and Illuminance	•	•				
Scatter and Leakage	•	•			•	•
DAP	•	•			•	•
Sensor Compatibility						
Accu-Gold Multi-Sensors (AGMS)	•	•	•	•		
Accu-Gold Ion Chambers (10x6)	•	•			•	•
DAP 60 Chamber	•	•			•	•
Accu-Gold Dose Diodes (DDX6)	•	•	•	•	•	•
DAPChkp+	•	•				
mA and mAs Sensor (90M)	•	•	•	•		
Light Sensor (AGLS)	•	•				
Computer Features						
Accu-Gold Windows Application		•		•		•
Waveforms	•	•	•	•	•	•
Excel Reporting		•		•		•
iOS and Android Applications		•		•		•
USB Connectivity		•		•		•
Wireless Connectivity		•		•		•